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**PERFORMANCE**

MERSILENE® sutura elicita un minimal initial inflammatory reaction in tissue followed by a gradual encapsulation of the suture by fibrous connective tissue. Subcutaneous tissue implantation in rats shows that no significant change in the retention of tensile strength of the suture occurs during the entire evaluation period of 180 days.

CONTRA-INDICATIONS

MERSILENE® sutures are sterilized by irradiation. Do not re-sterilize! Do not use if package is opened or damaged! Discard opened, unused sutures!

WARNINGS/PRECAUTIONS/INTERACTIONS

Users should be familiar with surgical procedures and techniques involving non absorbable sutures before employing MERSILENE® suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. As with any foreign body, prolonged contact of the suture with soft tissues, such as those found in the urinary or biliary tract, may result in calculus formation. Acceptable surgical practice should be followed for the management of infected or contaminated wounds.

In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance and the experience of the surgeon. The use of additional throws is particularly appropriate when knotting any monofilament suture.

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. All needles except bicolor needles labelled with MRI are magnetizable and should therefore not be used in an active magnetic field. Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in 'Sharps' containers.

**DESCRIPTION**

MERSILENE® suture is a monofilament or braided, non absorbable sterile surgical suture composed of poly (ethylene terephthalate). The empirical molecular formula of the polymer is  $(C_6H_{10}O_4)_n$ .

MERSILENE® is available undyed (white) or dyed green with D & C Green No.6, Colour Index No. 61565.

MERSILENE® complies with the requirements of the European Pharmacopoeia for Sterile Poly (ethylene Terephthalate) Suture and the United States Pharmacopoeia for Non Absorbable Surgical Sutures (except for an occasional slight increase in some soft tissue to bone tissue).

**APPLICATION**

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

**Instructions for use**

**MERSILENE\***  
 (POLYESTER)  
 STERILE SYNTETIC NON ABSORBABLE  
**SUTURE**

STERILITY

MERSILENE® sutures are sterilized by irradiation. Do not re-sterilize! Do not use if package is opened or damaged! Discard opened, unused sutures!

STORAGE

Recommended storage conditions: Below 25°C, away from moisture and direct heat. Do not use after expiry date!

SYMBOLS USED ON LABELLING

= Do not reuse  
 = Use until Year & Month

**CE 0086** = Sterile, unless package is damaged or opened. Method of Sterilization – Irradiation

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